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**Chemical behavior of:
Platinum-, Alumina Oxide-, and Graphite-Crucibles**

Platinum

Critical for Platinum is:

- Halogens (Cl_2 , F_2 , Br_2), aqua regia
- Li_2CO_3 , already before CO_2 -release
- SiC at app. 1000°C
- FeCl_3
- Be composites (evaporate just above melting point)
- HCL with oxidizer, for example (chromic acid, manganat, Iron (III)-salts, salt smelter)
- Reducing atmosphere
- Pb, Zn, Sn, Ag, Au, Hg, Li, Na, K, Sb, Bi, Ni, Fe, Steel, As, Si (composition creation)
- P, B
- Se above 320°C (cool down sample immediately after measurement and remove probe to prevent Selenium evaporation)
- Metal oxides with reducing substances such as C, organic compositions or H_2
- Sulfur, keying the surface
- Alkali hydroxide, - carbonates, - sulfates, - cyanides and - rhodanides at elevated temperatures
- KHSO_4 at elevated temperatures
- Carbon black or carbon above 1000°C
- SiO_2 under reducing conditions
- HBr, KCN-dilution at elevated temperatures
- High temperature resistant oxides above 1000°C

Platinum is not resistant against:

- Mixtures out of KNO_3 and NaOH at 700°C sealed from air
- Mixtures out of KOH and K_2S at 700°C sealed from air
- LiCl at 600°C
- Na_2O_2 at 500°C sealed from air
- MgCl_2 , $\text{Ba}(\text{NO}_3)_2$ at 700°C
- HBr, HJ, H_2O_2 (30%) and HNO_3 at 100°C
- KCl (damaging at the decomposition products crated during the melting, melting point 768°C)

Limited durable against:

- KHF_2 , LiF, NaCl at 900°C
- Mixtures of NaOH and NaNO_3 at 700°C sealed from air

Aluminumoxid

Critical for Aluminumoxid is:

- N₂ in presence of Carbon: creation of AlN, for this reason measurements with Carbon at elevated temperatures under nitrogen atmosphere are critical
- F₂: creation of AlF₃ and O₂
- Cl₂: creation of AlCl₃ above 700°C
- Sulfur: no reaction with liquid sulfur, if Carbon is present during the gas phase, creation of sulfides at elevated temperatures
- H₂S: under heat exposure creation of Al₂S up to 3%
- C: Creation of Carbides and Al at app. 1400°C
- HF: quantitative reaction to AlF₃ and H₂O at elevated temperatures
- CuSO₄: Diffusion trough crucible surface at app. 1000°C
- Compositions containing Silica, such as MOSi₂: in Inertia lead to contamination of Al₂O₃ at app. 1200°C; under air reaction at the contact point
- Metal fluorides: attacked by smelter, creation of (AlF₆)₃- anions and cryolite related salts
- Glass: Glass smelter dissolves Al₂O₃
- Hydrogen sulfate of the Alkali- and Alkaliearthoxide
- HCl: no reaction until 600°C, reaction at elevated temperatures which is increased when Carbon is present
- B₂O₃ or Borax: Smelter dilutes Al₂O₃, creation of Aluminoborates and -borides
- Alkali- and Alkaliearthoxide and their salts with volatile anions: smelter creates Aluminates or double oxides, important for Hydroxides, Nitrides, Nitrates, Carbonates, Peroxides, etc
- CaC₂: under heat creation of Al₄C₃
- PbO: Reaction at 700°C, also important for Lead oxide and Lead salt in combination with anions of volatile acids
- UO₃: Reaction start at 450°C, analog PbO
- MeO: Me = Fe₂₊, Co₂₊, Ni₂₊, etc.
- Alkali- Earthalcaliferrites : smelter dissolves Al₂O₃
- Zirconium composites with melting points between 800 and 1200°C : slow and weak reaction
- Some metal alloys such as: Fe with 4% Al

Graphite

Critical for Graphite is:

- O₂, reaction above 400°C
- Melted metals such as: Fe, Co, Ni, Na
- N₂ reaction starting from 1700°C (Creation of small amounts of Cyanides)
- Oxides (probable reduction when directly contacting)
- Water vapor
- F₂, Br₂, at room temperature
- Sulfur
- Si at app. 1400°C (formation of SiC)
- Chromic acid (aqueous)
- Cl-SO₃H
- SiO₂ → formation of SiC via the intermediate product SiO (technical production of SiC at temperatures above 1800°C, for that the reaction of SiO₂ and C certainly starts earlier)
- NO, NO₂
- Sulfur acid H₂SO₄ concentrated at arr. 150°C
- Saltpeter acid HNO₃, diluted at arr. 90°C
- SO₃ at arr. 100°C
- Danger of explosion with HClO₄
- NaOCl at arr. 50°C